Diagnostická váha iGET HOME B18

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Informace o produktu:

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User Manual – EN

Product information:

- Power supply 3xAAA batteries
- Maximum load capacity 180kg
- Accuracy +/- 100g
- Works with "iGET HOME" app analyzes 18 body values.
Weight / BMI / Fat % volume / Muscle mass (kg) / Body fat (kg) / Body fat index / Overweight level / Ideal weight / Weight control / Body fat / Fat free weight / Body water / Bone mass / Protein content / BMR / Metabolic age / Body type / Body score
- Trend line curve for history tracking
- More users
- Automatic on and off
- Low battery and overload indicator

Before first use:

- 1) Open the battery compartment on the bottom.
- 2) Insert the batteries, making sure the polarity is correct.
- 3) Close the battery compartment again.

The scales are for household use only, they cannot be used for commercial purposes.

Note: If the display shows the message "Lo", replace the batteries.

iGET HOME app

Download and install the iGET HOME app from Google Play or the App Store. Use the keyword "iGET HOME". Or use the following QR codes:

Android: Google Play

iOS: App Store

Adding a device:

- 1) Install the "iGET HOME" app from Google Play or the App Store.
- 2) Turn on bluetooth on your device.
- 3) Launch the "iGET HOME" app and follow the app's instructions. Create or log in to your account.
- 4) If you are using the app for the first time, you need to create an object first.
- 5) Then switch the application to the Devices tab
- 6) Turn on the scale, the scale will show 0.0 kg on the display. The Devices to be added screen will automatically appear on the application screen. Click on Add to add the device.
- 7) You will appear on the page where the device will be gradually added. When the pipe in the blue circle appears, click Next.

For full instructions on how to use the app, please visit: <https://www.iget.eu>

Safety warning:

- Do not expose the scale to open flames.
- Do not expose the scale to magnetic radiation.
- Avoid extreme temperatures (-10 °C, +40 °C).
- Do not disassemble the balance.
- In the event of any change in the dimensions of the balance due to overheating of the battery, stop using the device until the battery has been replaced with a new one.
- Do not use the scale if you have a pacemaker
- Do not use the scale if you are a pregnant woman
- Do not use the scale on slippery floors
- Keep in mind that this scale is not a medical device and its functions are set up to best suit the general population. Accurate advice on nutrition and fitness will be given to you by a professional.

Damage requiring professional repair:

If the following occurs, contact an authorized service center or supplier:

- Moisture has entered the scale.
- The scale is mechanically damaged.
- The scale is overheating.

For households: The symbol (crossed-out bin) on the product or in the accompanying documentation means that used electrical or electronic products must not be disposed of with municipal waste. In order to dispose of the product correctly, please take it to the designated collection points where it will be accepted free of charge. By properly disposing of this product, you will help conserve valuable natural resources and help prevent potential negative environmental and human health impacts that could result from improper waste disposal. For further details, contact your local authority or nearest collection point. Fines may be imposed in accordance with the regulations if this type of waste is disposed of incorrectly. Information for users on the disposal of electrical and electronic equipment (business and corporate use). Information for users on the disposal of electrical and electronic equipment in other countries outside the European Union: The above symbol (crossed out basket) is only valid in European Union countries. For proper disposal of electrical and electronic equipment, please ask your authorities or equipment dealer for detailed information. The symbol of the crossed-out bin on the product, packaging or printed materials represents everything. For warranty repairs of equipment, please contact your dealer. For technical problems and questions, contact your dealer. Observe the rules for working with electrical equipment. The user is not authorized to disassemble the equipment or replace any part of it. There is a risk of electric shock if the covers are opened or removed. Incorrect assembly and reconnection of the equipment also exposes you to the risk of electric shock.

**Declaration of conformity:**

INTELEK.CZ s.r.o. hereby declares that the BODY B18 type of radio equipment complies with Directive 2014/53/EU. The full text of the EU Declaration of Conformity is available on this website www.iget.eu.

The warranty period for the products is 24 months unless otherwise stated.

RoHS:

This device complies with the requirements of the Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment (Government Regulation No. 481/2012 Coll., as amended by Government Regulation No. 391/2016 Coll.) and thus with the requirements of Directive (EC) 2014/53/EU of the European Parliament and of the Council on the restriction of the use of certain hazardous substances in electrical and electronic equipment. The RoHS declaration can be downloaded at www.iget.eu.

This equipment can be used in the following countries:

AT	BE	CY	CZ	DK	EE	FI	MT	NL	PL	PT	SK
SI	ES	DE	GR	HU	IE	IT	LV	LT	GB	IS	LI
NO	CH	BG	HR	RU	RO	FR	SE	LU	TR		

The purchaser must ensure that the serial number is stated in the warranty card, as well as in the delivery note and proof of purchase. The Buyer acknowledges that if the serial number of the goods is not indicated on the purchase document, the delivery note and the warranty card and therefore it is not possible to verify that the goods are from the Seller by comparing the serial number on the goods and on their protective packaging with the serial number on the warranty card, the delivery note and the purchase document, the Service Centre is entitled to automatically and without further delay reject the claim for such goods. The Buyer is obliged to inform the persons to whom the goods are subsequently sold of this. Claim the scales from your Seller, unless otherwise stated on the website www.iget.eu. The warranty card is on the last pages of this manual or on the card itself.

Manufacturer / exclusive importer of iGET products into the EU:

INTELEK.CZ s.r.o. Olivova 2096/4,

110 00 Prague - CZ

Website: <http://www.iget.eu>

Support: <http://www.iget.eu/helpdesk>

Copyright © 2022 INTELEK.CZ s.r.o. All rights reserved.

Warranty card

Product:

Model:

Serial / Batch number (S/N):

Date of sale:

Warranty until:

Sale document number:

Seal of the store:

Information for Customers about Warranty Terms and Conditions:

Warranty shall apply to defects of purchased product, found during its warranty period. When purchasing a particular product, please consider carefully what properties and functions you expect your product to have. The fact that purchased product shall not satisfy your later requirements shall not give any reason for complaint. Please, prior to the first use of this product read carefully its operating instructions and follow these instructions consistently. We recommend you to keep original product packing during the warranty period for possible transport of the product in such a way so as to prevent any product damage during its transport and handling.

For the purpose of due application of any complaint about the product we recommend you to ask Seller to fill-in and confirm this Warranty Card including identification of serial number of purchased product, which shall comply with the serial number specified both in the document proving product purchase and delivery note and purchased product.

In this connection Customer hereby takes into account that unless this Warranty Card specifies serial number of purchased product, i.e., comparison of such serial number identified in the product with the serial number identified in this Warranty Card and document of product purchase/or delivery note shall not allow to confirm that particular product has come from Seller, in case of any doubts about the origin of such particular product Seller shall be entitled to automatically and without any other acts to reject the complaint of such goods.

Warranty shall not apply to:

product damage suffered by unqualified installation of firmware or its modification;

product damage by static electricity;

product use for any other purposes than designed, described in attached operating instructions;

damage to warranty seal of the producer, attached to the product;

product damage inconsistent with the operating instructions, technical standards or safety regulations valid in the Czech Republic;

product damage resulting from normal wear and tear;

discharged batteries, decrease in battery capacity caused by its wear and tear;

product defects caused by the use in dusty, polluted or in any other way dissatisfactory environment pursuant to attached operating instructions.

Wear and tear of the product and its components caused by their common use (discharged batteries, decreased capacity of accumulator, wear and tear of mechanical or displaying components, wear and tear of supply cables, headphones etc.) shall be subject neither to defect warranty nor quality warranty in the moment of take-over.

Furthermore, complaints about product parameters, not included in the operating instructions or in any commercial-technical documentation of Seller, related to the product, shall not be subject to quality warranty in the moment of take-over.

In the event of any changes in the dimensions of the battery, discontinue use, and until there is a change the battery with a new one.

Proof of product purchase, or delivery note, shall constitute an integral part of this Warranty Card.

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Copyright © 2022 INTELEK.CZ s.r.o. Všechna práva vyhrazena.
Copyright © 2022 INTELEK.CZ s.r.o. Všetky práva vyhradené.
Copyright © 2022 INTELEK.CZ s.r.o. Wszelkie prawa zastrzeżone.